




Product Data Sheet

Article Description	miro-absofoam	Foam Dressing
	<p>miro-absofoam is a sterile, semi-permeable, non-adhesive wound dressing with high absorption capacity. miro-absofoam maintains a moist environment to enhance the natural healing process. To be used during the granulation phase and for secretion absorption.</p> <p>miro-absofoam is permeable to gas but not to liquids and bacteria. The absorbent layer transfers exudate away from the wound and captures it without leaks or meltdown.</p> <p>The foam layer supports a high evaporation rate and cushions the wound. The dressing does not adhere to the wound, thereby reducing any pain associated with dressing removal preventing trauma to newly formed tissue. miro-absofoam has to be fixed with an appropriate secondary dressing.</p>	
	<p>INDICATIONS:</p> <ul style="list-style-type: none"> • exsuding acute and chronic wounds • pressure sores in all stages • ulcers of different origin 	

wound type	exsudation	to be used as
		primary dressing secondary dressing

Measurements


Item-N°	Size	Info	SU
24003	approx. 10 x 10cm		10 pieces

Composition of Goods

1. An open pores foam layer that provides a means for the escape of moisture vapour and offers a degree of cushioning to the wound, absorbs wound exudate
2. A vapour permeable, water and bacteria impermeable outer polyurethane layer

Storage Instructions

Storage: Store away from direct sunlight at ambient temperature (between 15°C and 25°C) and humidity.
Do not store next to heating and/or entrance doors.
„First in - first out“

■ necrotic wound
 ■ fibrinous coated wound
 ■ granulating wound
 ■ epithelising wound


Stretchable (in %)	Sterilization	Washable	CE (on single product)	Custom-Tariff-N°
---	Ethylenoxide	---	CE 0297	3005.9031
Shelf Life	Classification	Annex (EU-guideline 93/42)		Rule
5 years	I/s	V		4