




## Product Data Sheet

Article Description	miro-film frame	Sterile Wound Dressing
 	<p>miro-film is a sterile, self-adhesive wound dressing. The dressing is transparent and allows visual assessment of the wound without dressing change. It is very flexible and suitable for the protection of sites exposed to rubbing and shear stress. Used as a secondary dressing within the moist wound healing programme.</p> <p>miro-film is permeable to gas exchange but impermeable to water and bacteria. It does not adhere to the wound, thereby reducing any pain associated with dressing removal preventing trauma to newly formed tissue. Especially user-friendly is the white frame which helps to position and apply the dressing properly over the wound.</p> <p><b>INDICATIONS:</b></p> <ul style="list-style-type: none"> <li>• superficial chronic wounds (pressure sores, ulcers)</li> <li>• superficial acute wounds (burns, abrasions, donor sites)</li> <li>• Fixation of primary wound dressings (i.e. miro-gel, miro-sorb)</li> </ul>	
		

wound type	exsudation	to be used as
		primary dressing secondary dressing

### Measurements




Item-N°	Size	Info	SU
24708	approx. 6 x 8cm		50 pcs.
24712	approx. 10 x 12,5cm		50 pcs.
24730	approx. 20 x 30 cm		10 pcs.

### Composition of Goods

Thin polyurethan layer with hypoallergenic adhesive

### Storage Instructions

**Storage:** Store away from direct sunlight at ambient temperature (between 15°C and 25°C) and humidity.  
Do not store next to heating and/or entrance doors.  
„First in - first out“

■ necrotic wound ■ fibrinous coated wound ■ granulating wound ■ epithelising wound  
 low exsudation  moderate exsudation  heavily exsudation

Stretchable (in %)	Sterilization	Washable	CE (on single product)	Custom-Tariff-N°
---	Ethylenoxide	---	CE 0297	3005.1000
Shelf Life	Classification	Annex (EU-guideline 93/42)		Rule
5 years	I/s	V		4