


Product Data Sheet

Miro brand	miro-stretch transparent	<i>non-sterile, transparent fixation film</i>
Description	 <p>Non-sterile, length- and cross elastic wound dressing. miro-stretch transparent is used for securing dressings post operatively, and particularly over contoured areas, such as the nape of a joint and/or over skin folds. Highly conformable due to elastic PU foil. Permeable to air and vapour. Skinfriendly polyacrylate adhesive. Water resistant – very convenient for the patient.</p> <p style="text-align: right;">10m length per roll</p>	

Measurements

Item-N°	Size			PE
17105	5cm x 10m			1 piece
17110	10cm x 10m			1 piece
17115	15cm x 10m			1 piece

Composition of Goods

Backing material: transparent PU foil
Adhesive: Polyacrylate adhesive, hypo-allergenic, latex free

Storage Requirements:

Storage: Do not store extremely warm (over 25°C) or cold (under 10°C)
Avoid storage next to radiators or outer gates.
„First in - first out“.



Stretchable (in %)	Sterilization	Washable	CE (on single product)	Custom-Tariff-N°
---	Ethyleneoxide	---	yes	3005.1000
Shelf Life	Classification	Annex (EU-guideline 93/42)		Rule
5 years	I	I and VII		---